



## Pediatric Scrub Cap Instructions

### A. Cutting instructions:

1. Cut an 8" x WOF strip from a 1yd piece of fabric. Fold strip in half, lengthwise. Place hat band and tie piece on folds, as shown in figure 1a. Cut band and tie as one piece.

2. Cut another strip of fabric 7" x WOF, and fold in half, lengthwise. Place crown piece on fold as shown in figure 2a. One 7" strip will accommodate 4 crown pieces.



3. Repeat step 1 twice more.

\*Note - One yard of fabric is sufficient to make 3 caps (with one extra crown). If you have a smaller piece, and don't have enough for both the band and the crown, try making the crown out of a coordinating piece of fabric!

### B. Sewing instructions:

1. Fold band/tie piece with right sides together, and stitch a  $\frac{1}{4}$ " seam from each tie end to the notch, back-stitching at the notches. (see figure 1b.)

